



## ApoL2 Polyclonal Antibody

Cat No: HR1AP2794

For research use only

### Overview

Product Name	ApoL2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GenelD?Human?	APOL2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	APOL2; Apolipoprotein L2; Apolipoprotein L-II; ApoL-II
Background	apolipoprotein L2(APOL2) Homo sapiens This gene is a member of the apolipoprotein L gene family. The encoded protein is found in the cytoplasm, where it may affect the movement of lipids or allow the binding of lipids to organelles. Two transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008],