

## MEK2 Mouse mAb Cat No: HR1AM2316

For research use only

## Overview

Product Name	MEK2 Mouse mAb
Source	Mouse
Applications	WB
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	45kDakDa
GeneID?Human?	5605
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MAP2K2; MKK2; PRKMK2; Dual specificity mitogen-activated protein kinase kinase 2; MAP kinase kinase 2; MAPKK 2; ERK activator kinase 2; MAPK/ERK kinase 2
Background	MEK1 and MEK2, also called MAPK or Erk kinases, are dual-specificity protein kinases that function in a mitogen activated protein kinase cascade controlling cell growth and differentiation.