



## ApoF Polyclonal Antibody

Cat No: HR1AP2792

For research use only

### Overview

Product Name	ApoF Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GenID?Human?	APOF
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	APOF; Apolipoprotein F; Apo-F; Lipid transfer inhibitor protein; LTIP
Background	apolipoprotein F(APOF) Homo sapiens The product of this gene is one of the minor apolipoproteins found in plasma. This protein forms complexes with lipoproteins and may be involved in transport and/or esterification of cholesterol. [provided by RefSeq, Jul 2008],