



## VAMP-4 Polyclonal Antibody

Cat No: HR1AP8964

For research use only

### Overview

Product Name	VAMP-4 Polyclonal Antibody
Source	Rabbit
Applications	IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GenID?Human?	VAMP4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	VAMP4; Vesicle-associated membrane protein 4; VAMP-4
Background	vesicle associated membrane protein 4(VAMP4) Homo sapiens Synaptobrevins/VAMPs, syntaxins, and the 25-kD synaptosomal-associated protein SNAP25 are the main components of a protein complex involved in the docking and/or fusion of synaptic vesicles with the presynaptic membrane. The protein encoded by this gene is a member of the vesicle-associated membrane protein (VAMP)/synaptobrevin family. This protein may play a role in trans-Golgi network-to-endosome transport. [provided by RefSeq, Jul 2008].