



TAUT Polyclonal Antibody

Cat No: HR1AP8330

For research use only

Overview

Product Name	TAUT Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	70kDa
GenID?Human?	SLC6A6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SLC6A6; Sodium- and chloride-dependent taurine transporter; Solute carrier family 6 member 6
Background	solute carrier family 6 member 6(SLC6A6) <i>Homo sapiens</i> This gene encodes a multi-pass membrane protein that is a member of a family of sodium and chloride-ion dependent transporters. The encoded protein transports taurine and beta-alanine. There is a pseudogene for this gene on chromosome 21. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2013],