

Rabaptin-5? Polyclonal Antibody

Cat No: HR1AP4407

For research use only

Overview

Product Name	Rabaptin-5? Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	64kDa
GeneID?Human?	RABEP2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RABEP2; RABPT5B; Rab GTPase-binding effector protein 2; Rabaptin-5beta
Background	function:Plays a role in membrane trafficking and in homotypic early endosome fusion.,similarity:Belongs to the rabaptin family.,subunit:Heterodimer with RABGEF1. The dimer binds RAB5A that has been activated by GTP-binding.,