



CD6 Polyclonal Antibody

Cat No: HR1AP9106

For research use only

Overview

Product Name	CD6 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	73kDa
GeneID?Human?	CD6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CD6; T-cell differentiation antigen CD6; T12; TP120; CD antigen CD6
Background	CD6 molecule(CD6) Homo sapiens This gene encodes a protein found on the outer membrane of T-lymphocytes as well as some other immune cells. The encoded protein contains three scavenger receptor cysteine-rich (SRCR) domains and a binding site for an activated leukocyte cell adhesion molecule. The gene product is important for continuation of T cell activation. This gene may be associated with susceptibility to multiple sclerosis (PMID: 19525953, 21849685). Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2011],