



SEC22B Polyclonal Antibody

Cat No: HR1AP9473

For research use only

Overview

Product Name	SEC22B Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	23kDa
GenID?Human?	SEC22B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SEC22B; SEC22L1; Vesicle-trafficking protein SEC22b; ER-Golgi SNARE of 24 kDa; ERS-24; ERS24; SEC22 vesicle-trafficking protein homolog B; SEC22 vesicle-trafficking protein-like 1
Background	SEC22 homolog B, vesicle trafficking protein (gene/pseudogene)(SEC22B) Homo sapiens The protein encoded by this gene is a member of the SEC22 family of vesicle trafficking proteins. It seems to complex with SNARE and it is thought to play a role in the ER-Golgi protein trafficking. This protein has strong similarity to <i>Mus musculus</i> and <i>Cricetulus griseus</i> proteins.[provided by RefSeq, Sep 2009],