

NTF2 Polyclonal Antibody

Cat No: HR1AP5551

For research use only

Overview

Product Name	NTF2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	15kDa
GeneID?Human?	NUTF2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	NUTF2; NTF2; Nuclear transport factor 2; NTF-2; Placental protein 15; PP15
Background	nuclear transport factor 2(NUTF2) Homo sapiens This gene encodes a cytosolic factor that facilitates protein transport into the nucleus. The encoded protein is required for nuclear import of the small Ras-like GTPase, Ran which is involved in numerous cellular processes. This protein also interacts with the nuclear pore complex glycoprotein p62. [provided by RefSeq, Apr 2016],