



## NTF2 Polyclonal Antibody

Cat No: HR1AP5551

For research use only

### Overview

|                       |   |
|-----------------------|---|
| Product Name          | NTF2 Polyclonal Antibody  |
| Source                | Rabbit  |
| Applications          | WB,IHC-p,ELISA  |
| Species Reactivity    | Human,Mouse,Rat   |
| Recommended Dilutions |   |
| Immunogen             |   |
| Species               | Rabbit  |
| Storage               | -20°C/1 year  |
| Isotype               |   |
| Clonality             |   |
| Concentration         | 1 mg/ml   |
| Observed band         | 15kDa   |
| GenelD?Human?         | NUTF2   |
| Human Swiss-Prot No.  |   |
| Cellular localization |   |
| Alternative Names     | NUTF2; NTF2; Nuclear transport factor 2; NTF-2; Placental protein 15; PP15  |
| Background            | nuclear transport factor 2(NUTF2) Homo sapiens This gene encodes a cytosolic factor that facilitates protein transport into the nucleus. The encoded protein is required for nuclear import of the small Ras-like GTPase, Ran which is involved in numerous cellular processes. This protein also interacts with the nuclear pore complex glycoprotein p62. [provided by RefSeq, Apr 2016], |