



Amphiphysin I Polyclonal Antibody

Cat No: HR1AP2758

For research use only

Overview

Product Name	Amphiphysin I Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	76kDa
GenID?Human?	AMPH
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	AMPH; AMPH1; Amphiphysin
Background	amphiphysin(AMPH) Homo sapiens This gene encodes a protein associated with the cytoplasmic surface of synaptic vesicles. A subset of patients with stiff-man syndrome who were also affected by breast cancer are positive for autoantibodies against this protein. Alternate splicing of this gene results in two transcript variants encoding different isoforms. Additional splice variants have been described, but their full length sequences have not been determined. A pseudogene of this gene is found on chromosome 11.[provided by RefSeq, Nov 2010],