

Sec23B Polyclonal Antibody

Cat No: HR1AP5594

For research use only

Overview

Product Name	Sec23B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	86kDa
GeneID?Human?	SEC23B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SEC23B; Protein transport protein Sec23B; SEC23-related protein B
Background	Sec23 homolog B, coat complex II component(SEC23B) Homo sapiens The protein encoded by this gene is a member of the SEC23 subfamily of the SEC23/SEC24 family, which is involved in vesicle trafficking. The encoded protein has similarity to yeast Sec23p component of COPII. COPII is the coat protein complex responsible for vesicle budding from the ER. The function of this gene product has been implicated in cargo selection and concentration. Multiple alternatively spliced transcript variants have been identified in this gene. [provided by RefSeq, Feb 2010],