

CD28 (phospho Tyr218) Polyclonal Antibody

Cat No: HR1AP9143

For research use only

Overview

Product Name	CD28 (phospho Tyr218) Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GeneID?Human?	CD28
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CD28; T-cell-specific surface glycoprotein CD28; TP44; CD antigen CD28
Background	CD28 molecule(CD28) Homo sapiens The protein encoded by this gene is essential for T-cell proliferation and survival, cytokine production, and T-helper type-2 development. Several alternatively spliced transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Jul 2011],