

S1-Tag Mouse mAb

Cat No: HR1AM2306

For research use only

Overview

Product Name	S1-Tag Mouse mAb
Source	Mouse
Applications	WB
Species Reactivity	ALL
Recommended Dilutions	WB 1:5,000-10,000
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	N/AkDa
GeneID?Human?	N/A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	N/A
Background	S1 tag is an epitope tag composed of a nine-residue peptide, NANNPDWDF, derived from the hepatitis B virus preS1 region.