

# EPG5 Polyclonal Antibody

Cat No: HR1AP9311

For research use only

## Overview

Product Name	EPG5 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	
Observed band	170-292kDa
GeneID?Human?	EPG5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	EPG5; KIAA1632; Ectopic P granules protein 5 homolog
Background	ectopic P-granules autophagy protein 5 homolog(EPG5) Homo sapiens This gene encodes a large coiled coil domain-containing protein that functions in autophagy during starvation conditions. Mutations in this gene cause Vici syndrome. [provided by RefSeq, Aug 2015],