

# IKK? Polyclonal Antibody

Cat No: HR1AP3716

For research use only

## Overview

|                       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | IKK? Polyclonal Antibody                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Source                | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Applications          | WB,IHC-p,IF,ELISA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Species Reactivity    | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Immunogen             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Species               | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Isotype               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Clonality             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Observed band         | 43kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| GeneID?Human?         | IKBKG                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Cellular localization |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Alternative Names     | IKBKG; FIP3; NEMO; NF-kappa-B essential modulator; NEMO; FIP-3; Ikb kinase-associated protein 1; IKKAP1; Inhibitor of nuclear factor kappa-B kinase subunit gamma; I-kappa-B kinase subunit gamma; IKK-g                                                                                                                                                                                                                                                                                                                                                                                               |
| Background            | inhibitor of kappa light polypeptide gene enhancer in B-cells, kinase gamma(IKBKG) Homo sapiens This gene encodes the regulatory subunit of the inhibitor of kappaB kinase (IKK) complex, which activates NF-kappaB resulting in activation of genes involved in inflammation, immunity, cell survival, and other pathways. Mutations in this gene result in incontinentia pigmenti, hypohidrotic ectodermal dysplasia, and several other types of immunodeficiencies. A pseudogene highly similar to this locus is located in an adjacent region of the X chromosome. [provided by RefSeq, Mar 2016], |