

c-Fos Mouse mAb

Cat No: HR1AM2304

For research use only

Overview

Product Name	c-Fos Mouse mAb
Source	Mouse
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000-2,000 IHC 1:100-200
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	62kDa
GeneID?Human?	2353
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	AP1 antibody, Cellular oncogene c fos antibody, p55 antibody, proto oncogene c fos antibody
Background	c-Fos is a proto-oncogene that is the human homolog of the retroviral oncogene v-fos. c-fos is 380 amino acid protein with a basic leucine zipper region for dimerisation and DNA-binding and a transactivation domain at C-terminus.