

Wee 1 (phospho Ser53) Polyclonal Antibody

Cat No: HR1AP8627

For research use only

Overview

Product Name	Wee 1 (phospho Ser53) Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	WEE1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	WEE1; Wee1-like protein kinase; WEE1hu; Wee1A kinase
Background	WEE1 G2 checkpoint kinase(WEE1) Homo sapiens This gene encodes a nuclear protein, which is a tyrosine kinase belonging to the Ser/Thr family of protein kinases. This protein catalyzes the inhibitory tyrosine phosphorylation of CDC2/cyclin B kinase, and appears to coordinate the transition between DNA replication and mitosis by protecting the nucleus from cytoplasmically activated CDC2 kinase. [provided by RefSeq, Jul 2008],