

# Rsk-1 (phospho Thr573) Polyclonal Antibody

Cat No: HR1AP8193

For research use only

## Overview

|                       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | Rsk-1 (phospho Thr573) Polyclonal Antibody                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Source                | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Applications          | WB,IHC-p,IF,ELISA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Species Reactivity    | Human,Mouse,Rat                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Immunogen             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Species               | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Isotype               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Clonality             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Observed band         | 95kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| GeneID?Human?         | RPS6KA1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Cellular localization |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Alternative Names     | RPS6KA1; MAPKAPK1A; RSK1; Ribosomal protein S6 kinase alpha-1; S6K-alpha-1; 90 kDa ribosomal protein S6 kinase 1; p90-RSK 1; p90RSK1; p90S6K; MAP kinase-activated protein kinase 1a; MAPK-activated pro                                                                                                                                                                                                                                                                                                                                                                |
| Background            | ribosomal protein S6 kinase A1(RPS6KA1) Homo sapiens This gene encodes a member of the RSK (ribosomal S6 kinase) family of serine/threonine kinases. This kinase contains 2 nonidentical kinase catalytic domains and phosphorylates various substrates, including members of the mitogen-activated kinase (MAPK) signalling pathway. The activity of this protein has been implicated in controlling cell growth and differentiation. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008], |