

# HMG-1 Polyclonal Antibody

Cat No: HR1AP5343

For research use only

## Overview

|                       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|-----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | HMG-1 Polyclonal Antibody                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Source                | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Applications          | IF, WB, IHC-p, ELISA                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Species Reactivity    | Human, Mouse, Rat                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Immunogen             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Species               | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Isotype               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Clonality             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Observed band         | 25kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| GeneID?Human?         | HMGB1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Cellular localization |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Alternative Names     | HMGB1; HMG1; High mobility group protein B1; High mobility group protein 1; HMG-1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Background            | <p>high mobility group box 1(HMGB1) Homo sapiens This gene encodes a protein that belongs to the High Mobility Group-box superfamily. The encoded non-histone, nuclear DNA-binding protein regulates transcription, and is involved in organization of DNA. This protein plays a role in several cellular processes, including inflammation, cell differentiation and tumor cell migration. Multiple pseudogenes of this gene have been identified. Alternative splicing results in multiple transcript variants that encode the same protein. [provided by RefSeq, Sep 2015],</p> |