

CREB-1 (phospho Ser133) Polyclonal Antibody

Cat No: HR1AP2403

For research use only

Overview

Product Name	CREB-1 (phospho Ser133) Polyclonal Antibody
Source	Rabbit
Applications	IF, WB, IHC-p, IP, ELISA
Species Reactivity	Human, Mouse, Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	CREB1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CREB1; Cyclic AMP-responsive element-binding protein 1; CREB-1; cAMP-responsive element-binding protein 1
Background	cAMP responsive element binding protein 1(CREB1) Homo sapiens This gene encodes a transcription factor that is a member of the leucine zipper family of DNA binding proteins. This protein binds as a homodimer to the cAMP-responsive element, an octameric palindrome. The protein is phosphorylated by several protein kinases, and induces transcription of genes in response to hormonal stimulation of the cAMP pathway. Alternate splicing of this gene results in several transcript variants encoding different isoforms. [provided by RefSeq, Mar 2016],