

# Pragmin (phospho Tyr413) Polyclonal Antibody

Cat No: HR1AP6053

For research use only

## Overview

Product Name	Pragmin (phospho Tyr413) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	175kDa
GeneID?Human?	SGK223
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	SGK223; Tyrosine-protein kinase Sgk223; Sugen kinase 223
Background	homolog of rat pragma of Rnd2(SGK223) Homo sapiens This gene encodes an enzyme that belongs to the tyrosine protein kinase family. A similar protein in rat binds to Rho family GTPase 2 (Rnd2) and regulates neurite outgrowth via activation of Ras homolog gene family, member A (RhoA). [provided by RefSeq, Mar 2014],