



YAP (phospho Ser127) Polyclonal Antibody

Cat No: HR1AP2600

For research use only

Overview

Product Name	YAP (phospho Ser127) Polyclonal Antibody
Source	Rabbit
Applications	IF,WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	65kDa
GenelD?Human?	YAP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	YAP1; YAP65; Yorkie homolog; 65 kDa Yes-associated protein; YAP65
Background	Yes associated protein 1(YAP1) Homo sapiens This gene encodes a downstream nuclear effector of the Hippo signaling pathway which is involved in development, growth, repair, and homeostasis. This gene is known to play a role in the development and progression of multiple cancers as a transcriptional regulator of this signaling pathway and may function as a potential target for cancer treatment. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Aug 2013],