

# ZNF668 Polyclonal Antibody

Cat No: HR1AP8703

For research use only

## Overview

|                       |                                                                                                                                                                                  |
|-----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | ZNF668 Polyclonal Antibody                                                                                                                                                       |
| Source                | Rabbit                                                                                                                                                                           |
| Applications          | WB,IHC-p,IF,ELISA                                                                                                                                                                |
| Species Reactivity    | Human,Mouse                                                                                                                                                                      |
| Recommended Dilutions |                                                                                                                                                                                  |
| Immunogen             |                                                                                                                                                                                  |
| Species               | Rabbit                                                                                                                                                                           |
| Storage               | -20°C/1 year                                                                                                                                                                     |
| Isotype               |                                                                                                                                                                                  |
| Clonality             |                                                                                                                                                                                  |
| Concentration         | 1 mg/ml                                                                                                                                                                          |
| Observed band         | 68kDa                                                                                                                                                                            |
| GeneID?Human?         | ZNF668                                                                                                                                                                           |
| Human Swiss-Prot No.  |                                                                                                                                                                                  |
| Cellular localization |                                                                                                                                                                                  |
| Alternative Names     | ZNF668; Zinc finger protein 668                                                                                                                                                  |
| Background            | function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 16 C2H2-type zinc fingers., |