

ZNF541 Polyclonal Antibody

Cat No: HR1AP8855

For research use only

Overview

Product Name	ZNF541 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	83kDa
GeneID?Human?	ZNF541
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ZNF541; Zinc finger protein 541
Background	function:Component of some chromatin remodeling multiprotein complex that plays a role during spermatogenesis.,similarity:Contains 1 ELM2 domain.,similarity:Contains 1 SANT domain.,similarity:Contains 5 C2H2-type zinc fingers.,subunit:Interacts with HDAC1 and HSPA2.,