




ZNF600 Polyclonal Antibody

Cat No: HR1AP6082

For research use only

Overview

Product Name	ZNF600 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GenID?Human?	ZNF600
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ZNF600; Zinc finger protein 600
Background	function:May be involved in transcriptional regulation.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 20 C2H2-type zinc fingers.,



bioelsa