



ZNF592 Polyclonal Antibody

Cat No: HR1AP9190

For research use only

Overview

Product Name	ZNF592 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	160kDa
GenID?Human?	ZNF592
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ZNF592; KIAA0211; Zinc finger protein 592
Background	zinc finger protein 592(ZNF592) <i>Homo sapiens</i> This gene is thought to play a role in a complex developmental pathway and the regulation of genes involved in cerebellar development. Mutations in this gene have been associated with autosomal recessive spinocerebellar ataxia. [provided by RefSeq, Jan 2011],