

## Fyn Polyclonal Antibody

Cat No: HR1AP3493

For research use only

### Overview

|                       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|-----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | Fyn Polyclonal Antibody                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Source                | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Applications          | WB,IHC-p,ELISA                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Species Reactivity    | Human,Mouse,Rat                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Immunogen             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Species               | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Isotype               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Clonality             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Observed band         | 60kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| GeneID?Human?         | FYN                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Cellular localization |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Alternative Names     | FYN; Tyrosine-protein kinase Fyn; Proto-oncogene Syn; Proto-oncogene c-Fyn; Src-like kinase; SLK; p59-Fyn                                                                                                                                                                                                                                                                                                                                                                            |
| Background            | FYN proto-oncogene, Src family tyrosine kinase(FYN) Homo sapiens This gene is a member of the protein-tyrosine kinase oncogene family. It encodes a membrane-associated tyrosine kinase that has been implicated in the control of cell growth. The protein associates with the p85 subunit of phosphatidylinositol 3-kinase and interacts with the fyn-binding protein. Alternatively spliced transcript variants encoding distinct isoforms exist. [provided by RefSeq, Jul 2008], |