

ZNF262 Polyclonal Antibody

Cat No: HR1AP9099

For research use only

Overview

Product Name	ZNF262 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	175kDa
GeneID?Human?	ZMYM4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ZMYM4; KIAA0425; ZNF262; Zinc finger MYM-type protein 4; Zinc finger protein 262
Background	miscellaneous:The 3'-UTR region of the mRNA encoding this protein contains a motif called CDIR (for cell death inhibiting RNA) that binds HNRPD/AUF1 and HSPB1/HSP27. It is able to inhibit interferon-gamma induced apoptosis.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 9 MYM-type zinc fingers.,tissue specificity:Expressed at higher level in heart, skeletal muscle, kidney and liver.,