



ZNF446 Polyclonal Antibody

Cat No: HR1AP7852

For research use only

Overview

Product Name	ZNF446 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	49kDa
GeneID?Human?	ZNF446
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ZNF446; ZKSCAN20; Zinc finger protein 446; Zinc finger protein with KRAB and SCAN domains 20
Background	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 1 SCAN box domain.,similarity:Contains 3 C2H2-type zinc fingers.,