

## Polyphenol Oxidase (PPO) Activity Assay Kit

Cat No: HR3BC1265

For research use only

## Overview

Detection Method	Colorimetric method
Storage	2-8?
Instrument	Spectrophotometer(410 nm)
Assay Time	60 min
Validity	6
Assay Type	Enzyme Activity
Sample Type	Plant tissue
Synonyms	PPO : LIUEISU
Instrument	Spectrophotometer(410 nm)
Detection Principle	Polyphenol oxidase (PPO) can catalyze phenolic compounds into quinone substances. The latter has specific absorption at 410 nm. The activity of PPO can be calculated indirectly by measuring the OD value at 410 nm.
Reagents	
Labware	Tubes, Micropipette, Vortex mixer, 100? Water bath
Size	100 Assays
Sensitivity	
Detection Range	
Recovery Rate	