

Recombinant Human CDK1/CyclinA2 (N-GST tag)

Cat No:HR2R1252

For research use only

Overview

Quantity	10 ?g
Gene Symbol	CDK1/CyclinA2
Gene ID	N/A
Accession	NM_001786, NM_001237
Alternative Name	CDK1: CDC2 CDC28A; MGC111195; DKFZp686L20222 CyclinA2: CCN1; CCNA
Species	Human
Source	Insect cells
Description	CDK1 or Cell Division Control protein 1 is essential for the completion of START, the controlling event in the cell cycle that is required to initiate mitosis. CDK1 is a catalytic subunit of a protein kinase complex, called the M-Phase Promoting Factor that induces entry into mitosis and is universal among eukaryotes . Phosphorylation of Bcl-2 in G2/M phase-arrested cells following photodynamic therapy with hypericin involves a CDK1-mediated signal and delays the onset of apoptosis. Therapeutic potential of CDK inhibitor NU2058 in androgen-independent prostate cancer has also been demonstrated .
Functions	The specific activity of CDK1/CyclinA2 was determined to be 178 nmol /min/mg as per activity assay protocol.
Formulation	50mM Tris-HCl, pH 7.5, 150mM NaCl, 0.25mM DTT, 0.1mM EGTA, 0.1mM EDTA, 0.1mM PMSF, 25% glycerol.
Solubility	N/A
Appearance	Liquid
Molecular Weight	59
Purity	70% - 90%
Concentration	
Shipping Condition	Dry Ice
Storage Condition	Store product at ?70?C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.