

## Recombinant Human CDK2/CyclinA2 (N-GST tag)

Cat No:HR2R1255

For research use only

## Overview

Quantity	10 ?g
Gene Symbol	CDK2/CyclinA2
Gene ID	N/A
Accession	NM_001798, NM_001237
Alternative Name	CDK2: p33(CDK2) CyclinA2: CCN1; CCNA
Species	Human
Source	Insect cells
Description	CDK2 is a member of the Cyclin-Dependent Kinase family that is ubiquitously expressed. CDK2 is a catalytic subunit of the cyclin-dependent protein kinase complex, whose activity is restricted to the G1-S phase, and essential for cell cycle G1/S phase transition. CDK2 associates with and is regulated by the regulatory subunits of the complex including cyclin A or E, CDK inhibitor p21Cip1 (CDKN1A) and p27Kip1 (CDKN1B) . CDK2 phosphorylates multiple cellular substrates including SMAD3 and FOXO1. Phosphorylation of FOXO1 leads to its inhibition .?
Functions	The specific activity of CDK2/CyclinA2 was determined to be 1050 nmol/min/mg as per activity assay protocol.
Formulation	50mM Tris-HCl, pH 7.5, 150mM NaCl, 0.25mM DTT, 0.1mM EGTA, 0.1mM EDTA, 0.1mM PMSF, 25% glycerol.
Solubility	N/A
Appearance	Liquid
Molecular Weight	58
Purity	70% - 90%
Concentration	
Shipping Condition	Dry Ice
Storage Condition	Store product at ?70?C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.