

Recombinant Human CDK5/p35 (N-GST tag)

Cat No:HR2R1262

For research use only

Overview

Quantity	10 ?g
Gene Symbol	CDK5/p35
Gene ID	N/A
Accession	NM_004935, NM_003885
Alternative Name	CDK5: PSSALRE p35: CDK5R1; CDK5P35, CDK5R, NCK5A, p35, p35nck5a
Species	Human
Source	Insect cells
Description	CDK5 is a member of the Cyclin-Dependent Kinase family that is most abundant in the mammalian brain. Active form of CDK5, which has also been called neuronal cdc2-like kinase, is a heterodimer of CDK5 and a 25 kDa protein which is derived proteolytically from a 35 kDa brain and neuron-specific protein and is essential for the kinase activity of CDK5 . CDK5 has emerged as a crucial regulator of neuronal migration in the developing central nervous system. CDK5 phosphorylates a diverse list of substrates, implicating it in the regulation of a range of cellular processes - from adhesion and motility, to synaptic plasticity and drug addiction .
Functions	The specific activity of CDK5/p35 was determined to be 180 nmol /min/mg as per activity assay protocol.
Formulation	50mM Tris-HCl, pH 7.5, 150mM NaCl, 10mM glutathione, 0.1mM EDTA, 0.25mM DTT, 0.1mM PMSF, 25% glycerol.
Solubility	N/A
Appearance	Liquid
Molecular Weight	59
Purity	70% - 90%
Concentration	
Shipping Condition	Dry Ice
Storage Condition	Store product at ?70?C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.