

Recombinant Human CXCL2

Cat No:HR2R1314

For research use only

Overview

Quantity	1.0 ?g
Gene Symbol	CXCL2
Gene ID	2920
Accession	P19875
Alternative Name	C-X-C motif chemokine 2, Growth-regulated protein beta, Gro-beta, Macrophage inflammatory protein 2-alpha, MIP2-alpha, CXCL2, GRO2, GROB, MIP2A, SCYB2 Recombinant Human Growth Regulated Beta Protein (CXCL2)
Species	Human
Source	E. coli
Description	Growth regulated oncogene-gamma belongs to the family of chemotytic cytokines called chemokines. It is identical with MGSA (melanoma growth stimulatory activity) and the new designation is CXCL3. This factor is known mainly because of its chemotactic activity. GRO expression is inducible by serum or PDGF and/or by a variety of inflammatory mediators, such as IL-1 and TNF, in monocytes, fibroblasts, melanocytes and epithelial cells. In certain tumor cell lines, GRO is expressed constitutively. Similar to other alpha chemokines, the three GRO proteins are potent neutrophil attractants and activators. In addition, these chemokines are also active toward basophils. All three GROs can bind with high affinity to the IL-8 receptor type B.
Functions	Determined by its ability to chemoattract 293 transfected CXCR2 cells using a concentration range of 10-100 ng/mL.
Formulation	Recombinant GRO-beta/CXCL2 was lyophilized from a 0.2 ?m filtered concentrated (1.0 mg/mL) solution in 20 mM Tris, pH 8.0.
Solubility	A quick spin of the vial followed by reconstitution in distilled water to a concentration not less than 0.1 mg/mL. This solution can then be diluted into other buffers.
Appearance	Lyophilized Powder
Molecular Weight	8
Purity	>95% as determined by SDS-PAGE
Concentration	<1.0 EU/?g of recombinant protein as determined by the LAL method.
Shipping Condition	Ambient Temperature

Storage Condition	The lyophilized protein is stable for at least one year from date of receipt at -70°C. Upon reconstitution, this cytokine can be stored in working aliquots at 2° - 8°C for one month, or at -20°C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles.
-------------------	---

