

## Recombinant Human EPHB4 (N-GST tag)

Cat No:HR2R1364

For research use only

## Overview

Quantity	10 ?g
Gene Symbol	EPHB4
Gene ID	N/A
Accession	NM_004444
Alternative Name	HTK, MYK1, TYRO11, Mdk2
Species	Human
Source	Insect cells
Description	?EPHB4 is a member of the Eph family of receptor tyrosine kinases which play a role in regulating cell adhesion and cell movement during embryonic development . Furthermore, because tissue disorganization and abnormal cell adhesion, movement, and survival characterize the more advanced stages of cancer, the inappropriate functioning of Eph receptor could play a causal role in malignancy . EPHB4 is uniquely expressed in the vascular endothelial and endocardial cells and animal knockout studies with EPHB4 reveal a phenotype that is similar to Ephrin B2 gene knockdown. Furthermore, EPHB4 has been shown to be the major essential interaction partner of Ephrin B2 in angiogenesis.
Functions	The specific activity of EPHB4 was determined to be 350 nmol /min/mg as per activity assay protocol.
Formulation	50mM Tris-HCl, pH 7.5, 150mM NaCl, 0.25mM DTT, 0.1mM EGTA, 0.1mM EDTA, 0.1mM PMSF, 25% glycerol.
Solubility	N/A
Appearance	Liquid
Molecular Weight	50
Purity	70% - 90%
Concentration	
Shipping Condition	Dry Ice
Storage Condition	Store product at ?70?C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.