

Recombinant Human EPHB3 (N-GST tag)

Cat No:HR2R1363

For research use only

Overview

Quantity	10 ?g
Gene Symbol	EPHB3
Gene ID	N/A
Accession	NM_004443
Alternative Name	ETK2; HEK2; Cek10; Mdk5; TYRO6
Species	Human
Source	Insect cells
Description	EPHB3 is a member of the Ephrin receptor family and is expressed during embryonic development in multiple regions of the central nervous system. In adult brain, EPHB3 is expressed in the cerebellum, raphe pallidus, hippocampus, entorhinal cortex, and both motor and sensory cortices . EPHB3 is involved in the maintenance of mature neuronal connections and/or re-arrangement of synaptic connections during late stages of development. EPHB3 plays a role in the regulation of cell adhesion and migration, and the catalytic activity of EPHB3 is required for inhibition of integrin-mediated cell adhesion .?
Functions	The specific activity of EPBH3 was determined to be 26 nmol /min/mg as per activity assay protocol.
Formulation	50mM Tris-HCl, pH 7.5, 150mM NaCl, 0.25mM DTT, 0.1mM EGTA, 0.1mM EDTA, 0.1mM PMSF, 25% glycerol.
Solubility	N/A
Appearance	Liquid
Molecular Weight	68
Purity	70% - 90%
Concentration	
Shipping Condition	Dry Ice
Storage Condition	Store product at ?70?C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.