

Recombinant Human JAK1 (N-GST tag)

Cat No:HR2R1647

For research use only

Overview

Quantity	10 ?g
Gene Symbol	JAK1
Gene ID	N/A
Accession	NM_002227
Alternative Name	JAK1A, JAK1B
Species	Human
Source	Insect cells
Description	JAK1 is a member of protein-tyrosine kinases (PTK) characterized by the presence of a second phosphotransferase-related domain immediately N-terminal to the PTK domain. JAK1 bears all the hallmarks of a protein kinase, although its structure differs significantly from that of the PTK and threonine/serine kinase family members. JAK1 is a large, widely expressed membrane-associated phosphoprotein that is involved in the interferon-alpha/beta and -gamma signal transduction pathways. JAK1 plays an essential and nonredundant role in promoting biologic responses induced by a select subset of cytokine receptors, including those in which JAK utilization is thought to be nonspecific.
Functions	The specific activity of JAK1was determined to be 22nmol /min/mg as per activity assay protocol.
Formulation	50mM Tris-HCl, pH 7.5, 150mM NaCl, 10mM glutathione, 0.1mM EDTA, 0.25mM DTT, 0.1mM PMSF, 25% glycerol.
Solubility	N/A
Appearance	Liquid
Molecular Weight	108
Purity	70% - 90%
Concentration	
Shipping Condition	Dry Ice
Storage Condition	Store product at ?70?C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.