

Recombinant Human MAPKAPK2 (N-GST tag)

Cat No:HR2R1689

For research use only

Overview

Quantity	10 ?g
Gene Symbol	MAPKAPK2
Gene ID	N/A
Accession	NM_032960
Alternative Name	(2 isoforms)
Species	Human
Source	Insect cells
Description	MAPKAPK2 (MAPKAP kinase2) is a Ser/Thr protein kinase which is regulated via direct phosphorylation by p38 MAP kinase. In conjunction with p38 MAP kinase, MAPKAPK2 is known to be involved in many cellular processes including stress and inflammatory responses, nuclear export, gene expression regulation and cell proliferation. Heat shock protein HSP27 has been shown to be one of the substrates of MAPKAPK2 in vivo.?
Functions	The specific activity of MAPKAPK2 was determined to be 254 nmol /min/mg as per activity assay protocol.
Formulation	50mM Tris-HCl, pH 7.5, 150mM NaCl, 0.25mM DTT, 0.1mM EGTA, 0.1mM EDTA, 0.1mM PMSF, 25% glycerol.
Solubility	N/A
Appearance	Liquid
Molecular Weight	68
Purity	70% - 90%
Concentration	
Shipping Condition	Dry Ice
Storage Condition	Store product at ?70?C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.