

Recombinant Human Maspin (SerpineB5)

Cat No:HR2R1697

For research use only

Overview

Quantity	20 ?g
Gene Symbol	SERPINB5
Gene ID	5268
Accession	P36952
Alternative Name	Maspin, Peptidase inhibitor 5
Species	Human
Source	E. coli
Description	Mammary serine protease inhibitor (Maspin) is a non-inhibitory tumor suppressor serpin that is encoded by the SERPINB5 gene. Playing an important role in embryonic development in cell adhesion, Maspin has the ability to suppress the growth, invasion, and metastasis of breast and lung tumors. However, Maspin does not exhibit any serine protease inhibitor activity since it does not undergo the conformational transition from stressed to relaxed found in typical serpins. It is expressed in normal breast and prostatic epithelial cells, but downregulated in most breast carcinoma.
Functions	The ED(50) was determined by the dose-dependent proliferation of MDA-MB-468 cells and was found to be 0.05ng/mL.
Formulation	N/A
Solubility	N/A
Appearance	Lyophilized Powder
Molecular Weight	42.3
Purity	>95% as determined by SDS-PAGE
Concentration	1.0 EU/?g of recombinant protein as determined by the LAL method
Shipping Condition	Ambient Temperature
Storage Condition	The lyophilized protein is stable for at least one year from date of receipt at -70°C . Upon reconstitution, this cytokine can be stored in working aliquots at $2^{\circ} - 8^{\circ}\text{C}$ for one month, or at -20°C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles.