

Recombinant Human MUSK (N-GST tag)

Cat No:HR2R1771

For research use only

Overview

Quantity	10 ?g
Gene Symbol	MUSK
Gene ID	N/A
Accession	NM_005592
Alternative Name	MGC126323, MGC126324
Species	Human
Source	Insect cells
Description	MUSK is a receptor tyrosine kinase necessary for neuromuscular junction formation . MUSK gene expression is highly regulated during neuromuscular junction and it is involved in intercellular communication present on the surface of cells that are activated by specific protein ligands. MUSK members play a key role in growth and differentiation of those cell types. Agrin signals through MUSK to cluster acetylcholine receptors on the postsynaptic membrane of the neuromuscular junction. DOK7, a MUSK-interacting cytoplasmic protein, is essential for MUSK activation in cultured myotubes. MUSK also plays a critical role in the development of normal blood vessels .
Functions	The specific activity of MUSK was determined to be 6.7 nmol /min/mg as per activity assay protocol.
Formulation	50mM Tris-HCl, pH 7.5, 150mM NaCl, 0.25mM DTT, 0.1mM EGTA, 0.1mM EDTA, 0.1mM PMSF, 25% glycerol.
Solubility	N/A
Appearance	Liquid
Molecular Weight	63
Purity	70% - 90%
Concentration	
Shipping Condition	Dry Ice
Storage Condition	Store product at ?70?C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.