

## Recombinant Human MYLK4 (N-GST tag)

Cat No:HR2R1774

For research use only

## Overview

Quantity	10 ?g
Gene Symbol	MYLK4
Gene ID	N/A
Accession	BC132833
Alternative Name	SGK085
Species	Human
Source	Insect cells
Description	MYLK4 or myosin light chain kinase family, member 4 is Serine/Threonine protein kinases in which MYLK4 is one of four myosin light chain kinases in a comprehensive analysis of protein kinases encoded by the human genome.
Functions	The specific activity of MYLK4 was determined to be5 nmol /min/mg as per activity assay protocol.
Formulation	50mM Tris-HCl, pH 7.5, 150mM NaCl, 10mM glutathione, 0.1mM EDTA, 0.25mM DTT, 0.1mM PMSF, 25% glycerol.
Solubility	N/A
Appearance	Liquid
Molecular Weight	70
Purity	70% - 90%
Concentration	
Shipping Condition	Dry Ice
Storage Condition	Store product at ?70?C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.