

Recombinant Human NUAK1 (N-GST tag)

Cat No:HR2R1807

For research use only

Overview

Quantity	10 ?g
Gene Symbol	NUAK1
Gene ID	N/A
Accession	NM_014840
Alternative Name	ARK5; KIAA0537
Species	Human
Source	Insect cells
Description	NUAK1 or SNF1/AMP kinase-related kinase (SNARK) is a member of the NUAK family of SNF1-like kinase 1 that is also known as AMPK-related protein kinase 5 (ARK5). ARK5 is a tumor cell survival factor that is activated by AKT and acts as an ATM kinase under conditions of nutrient starvation . NUAK1 is highly expressed in heart and brain, and at lower levels in skeletal muscle, kidney, ovary, placenta, lung, and liver. NUAK1 is involved in tolerance to glucose starvation and suppresses Fas-induced apoptosis by phosphorylation of CASP6, thus suppressing the activation of the caspase and the subsequent cleavage of CFLAR .
Functions	The specific activity of NUAK1 was determined to be50 nmol /min/mg as per activity assay protocol.
Formulation	50mM Tris-HCl, pH 7.5, 150mM NaCl, 10mM glutathione, 0.1mM EDTA, 0.25mM DTT, 0.1mM PMSF, 25% glycerol.
Solubility	N/A
Appearance	Liquid
Molecular Weight	122
Purity	70% - 90%
Concentration	
Shipping Condition	Dry Ice
Storage Condition	Store product at ?70?C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.