

## Recombinant Human PAK1 (N-GST tag)

Cat No:HR2R1821

For research use only

## Overview

Quantity	10 ?g
Gene Symbol	PAK1
Gene ID	N/A
Accession	NM_002576
Alternative Name	PAKalpha, MGC130000, MGC130001
Species	Human
Source	Insect cells
Description	?PAK1 is a member of the p21-activated kinases (PAKs) which have been implicated in the regulation of cell morphology, motility and transformation. These serine/threonine kinases are activated by and are effectors of small GTPases, CDC42 and RAC. PAK1 belongs to the Group I PAKs which also includes PAK2 and PAK3. PAK1 is a key regulator of the actin cytoskeleton, adhesion and cell motility. Inactive dimeric PAK1 is mainly cytosolic and interaction with the activators Cdc42-GTP and Rac1-GTP stimulates the kinase at the sites of cellular protrusions forming adhesions to the extracellular matrix.
Functions	The specific activity of PAK1 was determined to be 74 nmol /min/mg as per activity assay protocol.
Formulation	50mM Tris-HCl, pH 7.5, 150mM NaCl, 10mM glutathione, 0.1mM EDTA, 0.25mM DTT, 0.1mM PMSF, 25% glycerol.
Solubility	N/A
Appearance	Liquid
Molecular Weight	95
Purity	70% - 90%
Concentration	
Shipping Condition	Dry Ice
Storage Condition	Store product at ?70?C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.