

Recombinant Human PIM2 (N-GST tag)

Cat No:HR2R1863

For research use only

Overview

Quantity	10 ?g
Gene Symbol	PIM2
Gene ID	N/A
Accession	NM_006875
Alternative Name	N/A
Species	Human
Source	Insect cells
Description	PIM2 is a serine threonine kinase that is present in all tissues, being most abundant in hematopoietic tissues, spleen, thymus, and peripheral blood leukocytes, as well as in testis, small intestine, and colon. It is highly expressed in human leukemic, lymphoma and colorectal adenocarcinoma cell lines. This suggests a role for PIM2 in proliferating cells as well as during meiosis. Similar to PIM1, PIM2 also acts as a pro-survival kinase and BAD protein is a legitimate PIM2 substrate.?
Functions	The specific activity of PIM2 was determined to be 305 nmol /min/mg as per activity assay protocol.
Formulation	50mM Tris-HCl, pH 7.5, 150mM NaCl, 0.25mM DTT, 0.1mM EGTA, 0.1mM EDTA, 0.1mM PMSF, 25% glycerol.
Solubility	N/A
Appearance	Liquid
Molecular Weight	61
Purity	70% - 90%
Concentration	
Shipping Condition	Dry Ice
Storage Condition	Store product at ?70?C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.