

Recombinant Human PKC alpha (N-GST tag)

Cat No:HR2R1872

For research use only

Overview

Quantity	10 ?g
Gene Symbol	PKCa
Gene ID	N/A
Accession	NM_002737
Alternative Name	AAG6, PKCA, PRKCA, PRKACA, MGC129900, MGC129901
Species	Human
Source	Insect cells
Description	PKCa is a member of the protein kinase C (PKC) family of serine- and threonine-specific protein kinases that can be activated by calcium and the second messenger diacylglycerol. PKC-alpha has been reported to play roles in many different cellular processes, such as cell adhesion, cell transformation, cell cycle checkpoint, and cell volume control. PKCa has been assigned to the chromosome region 17q22-q23.2 and has been identified as a fundamental regulator of cardiac contractility and Ca(2+) handling in myocytes .?
Functions	The specific activity of PKCa was determined to be 2900nmol/min/mg as per activity assay protocol
Formulation	50mM Tris-HCl, pH 7.5, 150mM NaCl, 0.25mM DTT, 0.1mM EGTA, 0.1mM EDTA, 0.1mM PMSF, 25% glycerol.
Solubility	N/A
Appearance	Liquid
Molecular Weight	103
Purity	70% - 90%
Concentration	
Shipping Condition	Dry Ice
Storage Condition	Store product at ?70?C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.