

Recombinant Human PRKG1 (N-GST tag)

Cat No:HR2R1896

For research use only

Overview

Quantity	10 ?g
Gene Symbol	PRKG1
Gene ID	N/A
Accession	NM_006258
Alternative Name	PGK, CGKI, PRKG1B, PRKGR1B, FLJ36117, MGC71944, cGKI-BETA, cGKI-alpha, DKFZp686K042
Species	Human
Source	Insect cells
Description	PRKG1 is a homodimer, with each monomer containing a regulatory cGMP-binding domain and a catalytic domain. By Northern blot analysis PRKG1 was shown to be expressed at highest levels in bladder, uterus, adrenal gland, and fallopian tube. PRKG1 plays an important stimulatory role in platelet activation. Expression of recombinant PRKG1 in a reconstituted cell model enhanced von Willebrand factor-induced activation of the platelet integrin alpha-IIb/beta-3. Prkg1 knockout mice showed impaired platelet responses to VWF or low doses of thrombin and prolonged bleeding time. Human platelet aggregation induced by VWF or low-dose thrombin was inhibited by PRKG1 inhibitors but enhanced by cGMP.
Functions	The specific activity of PRKG1 was determined to be 780 nmol /min/mg as per activity assay protocol.
Formulation	50mM Tris-HCl, pH 7.5, 150mM NaCl, 10mM glutathione, 0.1mM EDTA, 0.25mM DTT, 0.1mM PMSF, 25% glycerol.
Solubility	N/A
Appearance	Liquid
Molecular Weight	100
Purity	70% - 90%
Concentration	
Shipping Condition	Dry Ice
Storage Condition	Store product at ?70?C. For optimal storage, aliquot target into smaller quantities after centrifugation and store at recommended temperature. For most favorable performance, avoid repeated handling and multiple freeze/thaw cycles.