

Recombinant Mouse EGF

Cat No:HR2R2102

For research use only

Overview

Quantity	500 ?g
Gene Symbol	EGF
Gene ID	13645
Accession	P01132
Alternative Name	Epidermal growth factor
Species	Mouse
Source	E. coli
Description	<p>Epidermal growth factor (EGF) is a small mitogenic protein that is thought to be involved in mechanisms such as normal cell growth, oncogenesis, and wound healing. This protein shows both strong sequential and functional homology with human type-alpha transforming growth factor (hTGF alpha), which is a competitor for EGF receptor sites. EGF receptors are expressed in almost all types of tissues . Parietal endoderm, mature skeletal muscles, and hematopoietic tissues do not express the receptor. The EGF receptor, designated HER1, is a 170 kDa transmembrane glycoprotein with a length of 1186 amino acids. It is identical with a previously described glycoprotein called SA-7 (species antigen 7). The extracellular receptor domain has a length of 621 amino acids, including 11 glycosylated asparagine residues and 51 cysteine residues. This domain contains the EGF binding site and also binds mammalian TGF-alpha.</p>
Functions	The ED50 as determined by a cell proliferation assay using 3T3 cells was found to be ? 0.2 ng/mL.
Formulation	Lyophilized from a 0.2 ?m filtered solution in PBS, pH 7.5
Solubility	A quick spin of the vial followed by reconstitution in distilled water to a concentration not less than 0.1 mg/mL. This solution can then be diluted into other buffers
Appearance	Lyophilized Powder
Molecular Weight	6.2
Purity	>95% as determined by SDS-PAGE
Concentration	<1.0 EU/?g of recombinant protein as determined by the LAL method
Shipping Condition	Ambient Temperature
Storage Condition	The lyophilized protein is stable for at least one year from date of receipt at -70?C. Upon reconstitution, this cytokine can be stored in working aliquots at 2? - 8?C for one month, or at -20?C for six months, with a carrier protein without detectable loss of activity. Avoid repeated freeze/thaw cycles